


Product Name: Recombinant Human PCDHB9 Protein
Catalog #: PHH6152



Summary

Name	PCDB9_HUMAN
Purity	Greater than 85% as determined by reducing SDS-PAGE
Endotoxin level	Please contact with the lab for this information
Tag	C-hFc & C-His
Accession #	Q9Y5E1
Host	HEK293
Species	Human
Molecular Mass	40 kDa
Buffer	Phosphate buffered saline
Shipping	The product is shipped on dry ice/polar packs. Upon receipt, store it immediately at the temperature listed below.
Stability&Storage	Store at $\leq -70^{\circ}\text{C}$, stable for 6 months after receipt. Store at $\leq -70^{\circ}\text{C}$, stable for 3 months under sterile conditions after opening. Please minimize freeze-thaw cycles.
Sequence	KISENTAEGT AFRLEAQQDP DEGHNSIQNY TISSNSFFHI KISGSDEGMI YPELVLDKAL DREEQEELSL TLTALDGGSP SRSQTSTIRI VVLVDVNDNVP

image



Background

Alternative Names Protocadherin beta-9 (PCDH-beta-9) (Protocadherin-3H)

Note

For Research Use Only , Not for Diagnostic Use