


Product Name: Recombinant Human PSG2 Protein
Catalog #: PHH6169



Summary

Name	PSG2_HUMAN
Purity	Greater than 85% as determined by reducing SDS-PAGE
Endotoxin level	Please contact with the lab for this information
Tag	C-His
Accession #	P11465
Host	HEK293
Species	Human
Molecular Mass	34 kDa
Buffer	Phosphate buffered saline
Shipping	The product is shipped on dry ice/polar packs. Upon receipt, store it immediately at the temperature listed below.
Stability&Storage	Store at $\leq -70^{\circ}\text{C}$, stable for 6 months after receipt. Store at $\leq -70^{\circ}\text{C}$, stable for 3 months under sterile conditions after opening. Please minimize freeze-thaw cycles.
Sequence	QVTIEAQPPK VSEGKDVLLL VHNLPQNLTG YIWKYQGIRD LYHYITSYVV DGQIIYGPA YSGRETAYSN ASLLIQNVTR EDAGSYTLHI IKRGDGRGV TGYFTFTLYL ETPKPSISS NLNPREAMET VILTCPETP DTSYQWWMNG QSLPMTHRFQ LSETNRTLFL FGVTKYTAGP YECEIRNSGS ASRSDPVTLN LLHGPDLPRI HPSYTNYSRG DNLYLSCFAN SNPPAQYSWT INGKFQQSGQ NLFIPQITTK HSGLYVCSVR NSATGEESS SLTVKVSAST RIGLLPLLP T

image



Background

Alternative Names	Pregnancy-specific beta-1-glycoprotein 2 (PS-beta-G-2) (PSBG-2) (Pregnancy-specific glycoprotein 2) (Pregnancy-specific beta-1 glycoprotein E) (PS-beta-E)
--------------------------	--

Product Name: Recombinant Human PSG2 Protein
Catalog #: PHH6169



Note

For Research Use Only , Not for Diagnostic Use