


Product Name: Recombinant Human ERVH48-1 protein
Catalog #: PHH6188



Summary

Name	SUPYN_HUMAN
Purity	Greater than 85% as determined by reducing SDS-PAGE
Endotoxin level	Please contact with the lab for this information
Tag	N-MBP & C-His
Accession #	M5A8F1
Host	HEK293
Species	Human
Molecular Mass	56 kDa
Buffer	Phosphate buffered saline
Shipping	The product is shipped on dry ice/polar packs. Upon receipt, store it immediately at the temperature listed below.
Stability&Storage	Store at $\leq -70^{\circ}\text{C}$, stable for 6 months after receipt. Store at $\leq -70^{\circ}\text{C}$, stable for 3 months under sterile conditions after opening. Please minimize freeze-thaw cycles.
Sequence	APPSCRECYQ SLHYRGEMQQ YFTYHTHIER SCYGNLIEEC VESGKSYYKV KNLGVCGRN GAICPRGKQW LCFTKIGQWG VNTQVLEDIK REQIIAKAKA SKPTTPPENR PRHFHSFIQK L

image



SDS



WB

Background

Alternative Names	Suppressyn (Endogenous retrovirus group 48 member 1) (NDUFV3 antisense RNA 1) (endogenous retrovirus group Fb member 1)
--------------------------	---

Note

For Research Use Only , Not for Diagnostic Use