



Product Name: Recombinant Human SLC38A4 Protein
Catalog #: PHH6189



Summary

Name	S38A4_HUMAN
Purity	Greater than 85% as determined by reducing SDS-PAGE
Endotoxin level	Please contact with the lab for this information
Tag	C-hFc & C-His
Accession #	Q969I6
Host	HEK293
Species	Human
Molecular Mass	37 kDa
Buffer	Phosphate buffered saline
Shipping	The product is shipped on dry ice/polar packs. Upon receipt, store it immediately at the temperature listed below.
Stability&Storage	Store at $\leq -70^{\circ}\text{C}$, stable for 6 months after receipt. Store at $\leq -70^{\circ}\text{C}$, stable for 3 months under sterile conditions after opening. Please minimize freeze-thaw cycles.
Sequence	IEPDDESSG ESAPDSYIGI GNSEKAAMSS QFANEDTESQ KFLTNGFLGK KKLADYADEH HPGTTSFGMS

image



Background

Alternative Names	Sodium-coupled neutral amino acid transporter 4 (Amino acid transporter A3) (Na(+)-coupled neutral amino acid transporter 4) (Solute carrier family 38 member 4) (System A amino acid transporter 3) (System N amino acid transporter 3)
--------------------------	--

Note

For Research Use Only , Not for Diagnostic Use

Product Name: Recombinant Human SLC38A4 Protein
Catalog #: PHH6189

