


Product Name: Recombinant Human SLCO5A1 Protein
Catalog #: PHH6191



Summary

Name	SO5A1_HUMAN
Purity	Greater than 85% as determined by reducing SDS-PAGE
Endotoxin level	Please contact with the lab for this information
Tag	C-hFc & C-His
Accession #	Q9H2Y9
Host	HEK293
Species	Human
Molecular Mass	42 kDa
Buffer	Phosphate buffered saline
Shipping	The product is shipped on dry ice/polar packs. Upon receipt, store it immediately at the temperature listed below.
Stability&Storage	Store at $\leq -70^{\circ}\text{C}$, stable for 6 months after receipt. Store at $\leq -70^{\circ}\text{C}$, stable for 3 months under sterile conditions after opening. Please minimize freeze-thaw cycles.
Sequence	LGGINIPYTT GPSLTMPHRN LTGSCNVNCG CKIHEYEPVC GSDGITYFNP CLAGCVNSGN LSTGIRNYTE CTCVQSRQVI TPPTVGQRSQ LRVVIVKTYL NENGYAVSGK CKRTCNTLIP

image



SDS



WB

Background

Alternative Names	Solute carrier organic anion transporter family member 5A1 (Organic anion transporter polypeptide-related protein 4) (OATP-RP4) (OATPRP4) (Solute carrier family 21 member 15)
--------------------------	--

Note

For Research Use Only , Not for Diagnostic Use

Product Name: Recombinant Human SLC05A1 Protein
Catalog #: PHH6191

