


Product Name: Recombinant Human SLCO3A1 Protein
Catalog #: PHH6192



Summary

Name	SO3A1_HUMAN
Purity	Greater than 85% as determined by reducing SDS-PAGE
Endotoxin level	Please contact with the lab for this information
Tag	C-hFc & C-His
Accession #	Q9UIG8
Host	HEK293
Species	Human
Molecular Mass	38 kDa
Buffer	Lyophilized from Phosphate buffered saline
Shipping	The product is shipped on dry ice/polar packs. Upon receipt, store it immediately at the temperature listed below.
Stability&Storage	Store at $\leq -70^{\circ}\text{C}$, stable for 12 months after receipt. Store at $\leq -70^{\circ}\text{C}$, stable for 3 months under sterile conditions after opening. Please minimize freeze-thaw cycles.
Sequence	LGCDTG PVAG VTPYGNSTA PGSALDPYSP CNNNCECQTD SFTPVC GADG ITYLSACFAG CNSTNLTGCA CLTTVPAENA TVVPGKCPSP

image



Background

Alternative Names	Solute carrier organic anion transporter family member 3A1 (OATP3A1) (Organic anion transporter polypeptide-related protein 3) (OATP-RP3) (OATPRP3) (Organic anion-transporting polypeptide D) (OATP-D) (PGE1 transporter) (Sodium-independent organic anion transporter D) (Solute carrier family 21 member 11)
--------------------------	--

Note

Product Name: Recombinant Human SLC03A1 Protein
Catalog #: PHH6192



For Research Use Only , Not for Diagnostic Use