


Product Name: Recombinant Human SECTM1 Protein
Catalog #: PHH6196



Summary

Name	SCTM1_HUMAN
Purity	Greater than 85% as determined by reducing SDS-PAGE
Endotoxin level	Please contact with the lab for this information
Tag	C-His
Accession #	Q8WVN6
Host	HEK293
Species	Human
Molecular Mass	14 kDa
Buffer	Phosphate buffered saline
Shipping	The product is shipped on dry ice/polar packs. Upon receipt, store it immediately at the temperature listed below.
Stability&Storage	Store at $\leq -70^{\circ}\text{C}$, stable for 6 months after receipt. Store at $\leq -70^{\circ}\text{C}$, stable for 3 months under sterile conditions after opening. Please minimize freeze-thaw cycles.
Sequence	EGWDSPICTE GVSVSWGGEN TVMSCNISNA FSHVNIK LRA HGQESAIFNE VAPGYFSRDG WQLQVQGGVA QLVIKGARDS HAGLYMWHLV GHQRNNRQVT LEVSGAEPQS

image



Background

Alternative Names Secreted and transmembrane protein 1 (Protein K-12)

Note

For Research Use Only , Not for Diagnostic Use