


Product Name: Recombinant Human TSPAN6 Protein
Catalog #: PHH6199



Summary

Name	TSN6_HUMAN
Purity	Greater than 85% as determined by reducing SDS-PAGE
Endotoxin level	Please contact with the lab for this information
Tag	C-hFc & C-His
Accession #	O43657
Host	HEK293
Species	Human
Molecular Mass	37 kDa
Buffer	Phosphate buffered saline
Shipping	The product is shipped on dry ice/polar packs. Upon receipt, store it immediately at the temperature listed below.
Stability&Storage	Store at $\leq -70^{\circ}\text{C}$, stable for 6 months after receipt. Store at $\leq -70^{\circ}\text{C}$, stable for 3 months under sterile conditions after opening. Please minimize freeze-thaw cycles.
Sequence	SFKNNYEKAL KQYNSTGDYR SHAVDKIQNT LHCCGVTDYR DWTDTNYYSE KGFPKSCCKL EDCTPQRDAD

image



Background

Alternative Names	Tetraspanin-6 (Tspan-6) (A15 homolog) (Putative NF-kappa-B-activating protein 321) (T245 protein) (Tetraspanin TM4-D) (Transmembrane 4 superfamily member 6)
--------------------------	--

Note

For Research Use Only , Not for Diagnostic Use