



Product Name: Recombinant Human RYK Protein
Catalog #: PHH6212



Summary

Name	RYK_HUMAN
Purity	Greater than 85% as determined by reducing SDS-PAGE
Endotoxin level	Please contact with the lab for this information
Tag	C-hFc & C-His
Accession #	P34925
Host	HEK293
Species	Human
Molecular Mass	50 kDa
Buffer	Phosphate buffered saline
Shipping	The product is shipped on dry ice/polar packs. Upon receipt, store it immediately at the temperature listed below.
Stability&Storage	Store at $\leq -70^{\circ}\text{C}$, stable for 6 months after receipt. Store at $\leq -70^{\circ}\text{C}$, stable for 3 months under sterile conditions after opening. Please minimize freeze-thaw cycles.
Sequence	LLLALLPLLP APGAAAAPAP RPPELQSASA GPSVSLYLSE DEVRRLLIGLD AELYYVRNDL ISHYALSFSL LVPSETNFLH FTWHAKSKVE YKLGQVDNV LAMDMPPQVNI SVQGEVPRTL SVFRVELSCT GKVDSEVMIL MQLNLTVNSS KNFTVLNFKR RKMICYKKLEE VKTSALDKNT SRTIYDPVHA

image



Background

Alternative Names Tyrosine-protein kinase RYK (EC 2.7.10.1)

Note

For Research Use Only , Not for Diagnostic Use

Product Name: Recombinant Human RYK Protein
Catalog #: PHH6212

