

Summary

Name	SO6A1_HUMAN
Purity	Greater than 85% as determined by reducing SDS-PAGE
Endotoxin level	Please contact with the lab for this information
Tag	C-hFc & C-His
Accession #	Q86UG4
Host	HEK293
Species	Human
Molecular Mass	40 kDa
Buffer	Phosphate buffered saline
Shipping	The product is shipped on dry ice/polar packs. Upon receipt, store it immediately at the temperature listed below.
Stability&Storage	Store at $\leq -70^{\circ}\text{C}$, stable for 6 months after receipt. Store at $\leq -70^{\circ}\text{C}$, stable for 3 months under sterile conditions after opening. Please minimize freeze-thaw cycles.
Sequence	QFAGINEDYD GTGKLGNLTA PCNEKCRCSS SIYSSICGRD DIEYFSPCFA GCTYSKAQNNQ KKMYYNCSCI KEGLITADAE GDFIDARPGK CDAKCYKLPL

image



Background

Alternative Names	Solute carrier organic anion transporter family member 6A1 (Cancer/testis antigen 48) (CT48) (Gonad-specific transporter) (GST) (Organic anion-transporting polypeptide 6A1) (Organic anion-transporting polypeptide I) (OATP-I) (Solute carrier family 21 member 19)
--------------------------	---

Note

For Research Use Only , Not for Diagnostic Use

Product Name: Recombinant Human SLC06A1 Protein
Catalog #: PHH6214

