


Product Name: Recombinant Human SAA2 Protein
Catalog #: PHH6245




Summary

Name	SAA2_HUMAN
Purity	Greater than 85% as determined by reducing SDS-PAGE
Endotoxin level	Please contact with the lab for this information
Tag	N-MBP & C-His
Accession #	P0DJ19
Host	HEK293
Species	Human
Molecular Mass	54 kDa
Buffer	Phosphate buffered saline
Shipping	The product is shipped on dry ice/polar packs. Upon receipt, store it immediately at the temperature listed below.
Stability&Storage	Store at $\leq -70^{\circ}\text{C}$, stable for 6 months after receipt. Store at $\leq -70^{\circ}\text{C}$, stable for 3 months under sterile conditions after opening. Please minimize freeze-thaw cycles.
Sequence	RSFFSFLGEA FDGARDMWRA YSDMREANYI GSDKYFHARG NYDAAKRGPG GAWAAEVISN ARENIQRLTG RGAEDSLADQ AANKWGRSGR DPNHFRPAGL PEKY

image



SDS



WB

Background

Alternative Names Serum amyloid A-2 protein (SAA2) [Cleaved into: Amyloid A2 protein (AA2)]

Note

For Research Use Only , Not for Diagnostic Use