


Product Name: Recombinant Human VIPR2 Protein
Catalog #: PHH6251



Summary

Name	VIPR2_HUMAN
Purity	Greater than 85% as determined by reducing SDS-PAGE
Endotoxin level	Please contact with the lab for this information
Tag	C-His
Accession #	P41587
Host	HEK293
Species	Human
Molecular Mass	11 kDa
Buffer	Phosphate buffered saline
Shipping	The product is shipped on dry ice/polar packs. Upon receipt, store it immediately at the temperature listed below.
Stability&Storage	Store at $\leq -70^{\circ}\text{C}$, stable for 6 months after receipt. Store at $\leq -70^{\circ}\text{C}$, stable for 3 months under sterile conditions after opening. Please minimize freeze-thaw cycles.
Sequence	IQEETKCAE LLRSQTEKHK ACSGVWDNIT CWRPANVGET VTVPCPKVFS NFYSKAGNIS KNCTSDGWSE TFPDFVDACG

image



Background

Alternative Names	Vasoactive intestinal polypeptide receptor 2 (VIP-R-2) (Helodermin-preferring VIP receptor) (Pituitary adenylate cyclase-activating polypeptide type III receptor) (PACAP type III receptor) (PACAP-R-3) (PACAP-R3) (VPAC2)
--------------------------	---

Note

For Research Use Only , Not for Diagnostic Use