



Product Name: Recombinant Human TMEFF1 Protein
Catalog #: PHH6254



Summary

Name	TEFF1_HUMAN
Purity	Greater than 85% as determined by reducing SDS-PAGE
Endotoxin level	Please contact with the lab for this information
Tag	C-His
Accession #	Q8IYR6
Host	HEK293
Species	Human
Molecular Mass	13 kDa
Buffer	Phosphate buffered saline
Shipping	The product is shipped on dry ice/polar packs. Upon receipt, store it immediately at the temperature listed below.
Stability&Storage	Store at $\leq -70^{\circ}\text{C}$, stable for 6 months after receipt. Store at $\leq -70^{\circ}\text{C}$, stable for 3 months under sterile conditions after opening. Please minimize freeze-thaw cycles.
Sequence	GPCKYKAECD EDAENVGCVC NIDCSGYSFN PVCASDGSSY NNPCFVREAS CIKQEIDIR HLGHCTDIDD TSLGKKDDG LQYRPDVKDA

image



Background

Alternative Names	Tomoregulin-1 (TR-1) (H7365) (Transmembrane protein with EGF-like and one follistatin-like domain)
--------------------------	--

Note

For Research Use Only , Not for Diagnostic Use