


Product Name: Recombinant Human VWC2L Protein
Catalog #: PHH6292



Summary

Name	VWC2L_HUMAN
Purity	Greater than 85% as determined by reducing SDS-PAGE
Endotoxin level	Please contact with the lab for this information
Tag	C-His
Accession #	B2RUY7
Host	HEK293
Species	Human
Molecular Mass	23 kDa
Buffer	Phosphate buffered saline
Shipping	The product is shipped on dry ice/polar packs. Upon receipt, store it immediately at the temperature listed below.
Stability&Storage	Store at $\leq -70^{\circ}\text{C}$, stable for 6 months after receipt. Store at $\leq -70^{\circ}\text{C}$, stable for 3 months under sterile conditions after opening. Please minimize freeze-thaw cycles.
Sequence	AAISHEDYPA DEGDQISSND NLIFDDYRGK GCVDDSGFVY KLGERFFPGH SNCPCVCALD GPVCDQPECP KIHPKCTKVE HNGCCPECKE VKNFCEYHGK NYKILEEFKP SPCEWCRCEP SNEVHCVVAD CAVPECVNPV YEPEQCCPVC KNGPNCFAGT TIIPAGIEVK VDECNICHCH NGDWWKPAQC SKRECQGKQT V

image



Background

Alternative Names von Willebrand factor C domain-containing protein 2-like (Brorin-like)

Note

For Research Use Only , Not for Diagnostic Use

Product Name: Recombinant Human VWC2L Protein
Catalog #: PHH6292

