



**Product Name: Recombinant Human ZDHHC14 Protein**  
**Catalog #: PHH6298**



## Summary

|                              |  |
|------------------------------|--|
| <b>Name</b>                  | ZDH14_HUMAN  |
| <b>Purity</b>                | Greater than 85% as determined by reducing SDS-PAGE  |
| <b>Endotoxin level</b>       | Please contact with the lab for this information   |
| <b>Tag</b>                   | C-His  |
| <b>Accession #</b>           | Q8IZN3   |
| <b>Host</b>                  | HEK293   |
| <b>Species</b>               | Human  |
| <b>Molecular Mass</b>        | 16 kDa   |
| <b>Buffer</b>                | Phosphate buffered saline  |
| <b>Shipping</b>              | The product is shipped on dry ice/polar packs. Upon receipt, store it immediately at the temperature listed below.   |
| <b>Stability&amp;Storage</b> | Store at $\leq -70^{\circ}\text{C}$ , stable for 6 months after receipt. Store at $\leq -70^{\circ}\text{C}$ , stable for 3 months under sterile conditions after opening. Please minimize freeze-thaw cycles. |
| <b>Sequence</b>              | QSQSDMCDQD QCIQSTKFVL QAAATPLLQS EPSLTSDELH LPGKPGLGTP<br>CASLTLGPPT PPASMPNLAE ATLADVMPRK DEHMGHQFLT PDEAPSPRL<br>LAAGSPLAHS RTMHVLGLAS QDSLHEDSVR  |

## image



## Background

|                          |   |
|--------------------------|---|
| <b>Alternative Names</b> | Palmitoyltransferase ZDHHC14 (EC 2.3.1.225) (DHHC domain-containing cysteine-rich protein 14) (DHHC-14) (NEW1 domain-containing protein) (NEW1CP) (Zinc finger DHHC domain-containing protein 14) |
|--------------------------|---|

## Note

For Research Use Only , Not for Diagnostic Use

**Product Name: Recombinant Human ZDHHC14 Protein**  
**Catalog #: PHH6298**

---

