



Product Name: Recombinant Human PSG3 Protein
Catalog #: PHH6320



Summary

Name	PSG3_HUMAN
Purity	Greater than 85% as determined by reducing SDS-PAGE
Endotoxin level	Please contact with the lab for this information
Tag	N-MBP & C-His
Accession #	Q16557
Host	HEK293
Species	Human
Molecular Mass	61 kDa
Buffer	Phosphate buffered saline
Shipping	The product is shipped on dry ice/polar packs. Upon receipt, store it immediately at the temperature listed below.
Stability&Storage	Store at $\leq -70^{\circ}\text{C}$, stable for 12 months after receipt. Store at $\leq -70^{\circ}\text{C}$, stable for 3 months under sterile conditions after opening. Please minimize freeze-thaw cycles.
Sequence	EPTKVSKGKD VLLLVHNLQP NLAGYIWYKG QMKDLYHYIT SYVVDGQIII YGPAYSGRET VYSNASLLIQ NVTREDAGSY TLHIVKRGDG TRGETGHFTF TLYLETPKPS ISSSNLYPRE DMEAVSLTCD PETPDASYLW WMINGQSLPMT HSLQLSKNKR

image



Background

Alternative Names Pregnancy-specific beta-1-glycoprotein 3 (PS-beta-G-3) (PSBG-3) (Pregnancy-specific glycoprotein 3) (Carcinoembryonic antigen SG5)

Note

For Research Use Only , Not for Diagnostic Use

Product Name: Recombinant Human PSG3 Protein
Catalog #: PHH6320

