


Product Name: Recombinant Human SVOP Protein
Catalog #: PHH6321



Summary

Name	SVOP_HUMAN
Purity	Greater than 85% as determined by reducing SDS-PAGE
Endotoxin level	Please contact with the lab for this information
Tag	C-hFc & C-His
Accession #	Q8N4V2
Host	HEK293
Species	Human
Molecular Mass	37 kDa
Buffer	Phosphate buffered saline
Shipping	The product is shipped on dry ice/polar packs. Upon receipt, store it immediately at the temperature listed below.
Stability&Storage	Store at $\leq -70^{\circ}\text{C}$, stable for 12 months after receipt. Store at $\leq -70^{\circ}\text{C}$, stable for 3 months under sterile conditions after opening. Please minimize freeze-thaw cycles.
Sequence	LPVVKFRRTG ESARSEDDTA SGEHEVQIEG VHVGLEAVEL DDGAAVPKEF ANPTDDTFMV EDAVEAIGFG

image



Background

Alternative Names Synaptic vesicle 2-related protein (SV2-related protein)

Note

For Research Use Only , Not for Diagnostic Use