


Product Name: Recombinant Human MUC21 Protein
Catalog #: PHH6343



Summary

Name	MUC21_HUMAN
Purity	Greater than 85% as determined by reducing SDS-PAGE
Endotoxin level	Please contact with the lab for this information
Tag	C-hFc & C-His
Accession #	Q5SSG8
Host	HEK293
Species	Human
Molecular Mass	43 kDa
Buffer	Phosphate buffered saline
Shipping	The product is shipped on dry ice/polar packs. Upon receipt, store it immediately at the temperature listed below.
Stability&Storage	Store at $\leq -70^{\circ}\text{C}$, stable for 6 months after receipt. Store at $\leq -70^{\circ}\text{C}$, stable for 3 months under sterile conditions after opening. Please minimize freeze-thaw cycles.
Sequence	ANTGSSVISS GASTATNSGS SVTSSGVSTA TISGSSVTSN GVSIVTNSEF HTTSSGISTA TNSEFSTVSS GISIATNSES STTSSGASTA TNSESSTPSS GASTATNSDS STTSSGASTA TNSDSSTTSS EASTATNSES STTSSGASTA TNSESSTVSS

image



Background

Alternative Names Mucin-21 (MUC-21) (Epiglycanin)

Note

For Research Use Only , Not for Diagnostic Use

Product Name: Recombinant Human MUC21 Protein
Catalog #: PHH6343

