


Product Name: Recombinant Human TNFRSF18 Protein
Catalog #: PHH6383



Summary

Name	TNR18_HUMAN
Purity	Greater than 85% as determined by reducing SDS-PAGE
Endotoxin level	Please contact with the lab for this information
Tag	C-His
Accession #	Q9Y5U5
Host	HEK293
Species	Human
Molecular Mass	16 kDa
Buffer	Phosphate buffered saline
Shipping	The product is shipped on dry ice/polar packs. Upon receipt, store it immediately at the temperature listed below.
Stability&Storage	Store at $\leq -70^{\circ}\text{C}$, stable for 6 months after receipt. Store at $\leq -70^{\circ}\text{C}$, stable for 3 months under sterile conditions after opening. Please minimize freeze-thaw cycles.
Sequence	GPGCGPGRLL LGTGTDARCC RVHTTRCCRD YPGEECCSEW DCMCVQPEFH CGDPCCTTCR HHPCPPGQGV QSQGKFSFGF QCIDCASGTF SGGHEGHCKP WTDCTQFGFL TVFPGNKTHN AVCVPGSPPA

image



Background

Alternative Names	Tumor necrosis factor receptor superfamily member 18 (Activation-inducible TNFR family receptor) (Glucocorticoid-induced TNFR-related protein) (CD antigen CD357)
--------------------------	---

Note

For Research Use Only , Not for Diagnostic Use

Product Name: Recombinant Human TNFRSF18 Protein
Catalog #: PHH6383

