



Product Name: Recombinant Human UNC5A Protein
Catalog #: PHH6386



Summary

Name	UNC5A_HUMAN
Purity	Greater than 85% as determined by reducing SDS-PAGE
Endotoxin level	Please contact with the lab for this information
Tag	C-His
Accession #	Q6ZN44
Host	HEK293
Species	Human
Molecular Mass	32 kDa
Buffer	Phosphate buffered saline
Shipping	The product is shipped on dry ice/polar packs. Upon receipt, store it immediately at the temperature listed below.
Stability&Storage	Store at $\leq -70^{\circ}\text{C}$, stable for 6 months after receipt. Store at $\leq -70^{\circ}\text{C}$, stable for 3 months under sterile conditions after opening. Please minimize freeze-thaw cycles.
Sequence	VANPVPGANP DLLPHFLVEP EDVYIVKNKP VLLVCKAVPA TQIFFKCNGE WVRQVDHVIE RSTDGSSGLP TMEVRINVSQ QQVEKVFGL EYWCQCVAWS SSGTTKSQKA YIRIAYLRKN FEQEPLAKEV SLEQGIVLPC RPPEGIPPAE VEWLRNEDLV DPSLDPNVYI TREHSLVVRQ ARLADTANYT CVAKNIVARR RSASAIVIVY VDGSWSPWSK WSACGLDCTH WRSRECSIPA PRNGGEECQG TDLDRNCTS DLCVHTASGP

image



Background

Alternative Names Netrin receptor UNC5A (Protein unc-5 homolog 1) (Protein unc-5 homolog A)

Product Name: Recombinant Human UNC5A Protein
Catalog #: PHH6386



Note

For Research Use Only , Not for Diagnostic Use