



Product Name: Recombinant Human TAC3 Protein
Catalog #: PHH6387



Summary

Name	TKNK_HUMAN
Purity	Greater than 85% as determined by reducing SDS-PAGE
Endotoxin level	Please contact with the lab for this information
Tag	N-9×His
Accession #	Q9UHF0
Host	HEK293
Species	Human
Molecular Mass	13 kDa
Buffer	Phosphate buffered saline
Shipping	The product is shipped on dry ice/polar packs. Upon receipt, store it immediately at the temperature listed below.
Stability&Storage	Store at $\leq -70^{\circ}\text{C}$, stable for 6 months after receipt. Store at $\leq -70^{\circ}\text{C}$, stable for 3 months under sterile conditions after opening. Please minimize freeze-thaw cycles.
Sequence	AVCKEPQEEV VPGGGRSKRD PDLYQLLQRL FKSHSSLEGL LKALSQASTD PKESTSPEKR DMHDFVGLM GKRSVQPDSP TDVNQENVPS FGILKYPPRA

image



Background

Alternative Names Tachykinin-3 (ZNEUROK1) [Cleaved into: Neurokinin-B (NKB) (Neuromedin-K)]

Note

For Research Use Only , Not for Diagnostic Use