


Product Name: Recombinant Human PSG5 Protein
Catalog #: PHH6431



Summary

Name	PSG5_HUMAN
Purity	Greater than 85% as determined by reducing SDS-PAGE
Endotoxin level	Please contact with the lab for this information
Tag	C-His
Accession #	Q15238
Host	HEK293
Species	Human
Molecular Mass	35 kDa
Buffer	Phosphate buffered saline
Shipping	The product is shipped on dry ice/polar packs. Upon receipt, store it immediately at the temperature listed below.
Stability&Storage	Store at $\leq -70^{\circ}\text{C}$, stable for 6 months after receipt. Store at $\leq -70^{\circ}\text{C}$, stable for 3 months under sterile conditions after opening. Please minimize freeze-thaw cycles.
Sequence	QVTIEALPPK VSEGKDVLLL VHNLPQNLG YIWYKQLMD LYHYITSYVV DGQINIYGPA YTGRETVYSN ASLLIQNVTR EDAGSYTLHI IKRGDRTRGV TGYFTFNLYL KLPKPYITIN NSKPRENKDV LAFTCEPKSE NYTYIWWLNG QSLPVSPRVK RPIENRILIL PSVTRNETGP YECEIRDRDG GMRSDPVTLN VLYGPDLPIS YPSFTYYRSG ENLYLSCFAE SNPPAEYFWT INGKFQQSGQ KLSIPQITTK HRGLYTCSVR NSATGKESK SMTVEVSAPS GIGRLPLLNP I

image



Background

Alternative Names	Pregnancy-specific beta-1-glycoprotein 5 (PS-beta-G-5) (PSBG-5) (Pregnancy-specific glycoprotein 5) (Fetal liver non-specific cross-reactive antigen 3) (FL-NCA-3)
--------------------------	--

Product Name: Recombinant Human PSG5 Protein
Catalog #: PHH6431



Note

For Research Use Only , Not for Diagnostic Use