


Product Name: Recombinant Human RECK Protein
Catalog #: PHH6443



Summary

Name	RECK_HUMAN
Purity	Greater than 85% as determined by reducing SDS-PAGE
Endotoxin level	Please contact with the lab for this information
Tag	N-MBP & C-His
Accession #	O95980
Host	HEK293
Species	Human
Molecular Mass	57 kDa
Buffer	Phosphate buffered saline
Shipping	The product is shipped on dry ice/polar packs. Upon receipt, store it immediately at the temperature listed below.
Stability&Storage	Store at $\leq -70^{\circ}\text{C}$, stable for 6 months after receipt. Store at $\leq -70^{\circ}\text{C}$, stable for 3 months under sterile conditions after opening. Please minimize freeze-thaw cycles.
Sequence	CMSKDPCNPN PCQKNQRCIP KPQVCLTTFD KFGCSQYECV PRQLACDQVQ DPVCDTDHME HNNLCTLYQR GKSLSYKGPC QPFCRATEPV CGHNGETYSS VCAAYSDRVA VDYYGDCQAV GVLSEHSSVA

image



Background

Alternative Names	Reversion-inducing cysteine-rich protein with Kazal motifs (hRECK) (Suppressor of tumorigenicity 15 protein)
--------------------------	---

Note

For Research Use Only , Not for Diagnostic Use