


Product Name: Recombinant Human PROM2 Protein
Catalog #: PHH6445



Summary

Name	PROM2_HUMAN
Purity	Greater than 85% as determined by reducing SDS-PAGE
Endotoxin level	Please contact with the lab for this information
Tag	C-hFc & C-His
Accession #	Q8N271
Host	HEK293
Species	Human
Molecular Mass	40 kDa
Buffer	Phosphate buffered saline
Shipping	The product is shipped on dry ice/polar packs. Upon receipt, store it immediately at the temperature listed below.
Stability&Storage	Store at $\leq -70^{\circ}\text{C}$, stable for 6 months after receipt. Store at $\leq -70^{\circ}\text{C}$, stable for 3 months under sterile conditions after opening. Please minimize freeze-thaw cycles.
Sequence	HLDINQYTNK LRQELQSLKV DTQSLDLLSS AARRDLEALQ SSSLQRIHYP DFLVQIQRPV VKTSMEQLAQ ELQGLAQDQD NSVLGQRLQE EAQGLRNLHQ EKVVPQQSLV

image



Background

Alternative Names Prominin-2 (PROM-2) (Prominin-like protein 2) (hPROML2)

Note

For Research Use Only , Not for Diagnostic Use