


Product Name: Recombinant Human ULBP3 Protein
Catalog #: PHH6508



Summary

Name	ULBP3_HUMAN
Purity	Greater than 85% as determined by reducing SDS-PAGE
Endotoxin level	Please contact with the lab for this information
Tag	C-His
Accession #	Q9BZM4
Host	HEK293
Species	Human
Molecular Mass	23 kDa
Buffer	Phosphate buffered saline
Shipping	The product is shipped on dry ice/polar packs. Upon receipt, store it immediately at the temperature listed below.
Stability&Storage	Store at $\leq -70^{\circ}\text{C}$, stable for 6 months after receipt. Store at $\leq -70^{\circ}\text{C}$, stable for 3 months under sterile conditions after opening. Please minimize freeze-thaw cycles.
Sequence	AHSLWYNFTI IHLPRHGQQW CEVQSQVDQK NFLSYDCGSD KVLSMGHLEE QLYATDAWGK QLEMLREVGQ RLRLELADTE LEDFTPSGPL TLQVRMSCEC EADGYIRGSW QFSFDGRKFL LFDSNNRKWT VVHAGARRMK EKWEKDSGLT TFFKMVSMRD CKSWLRDFLM HRKKRLEPTA

image



Background

Alternative Names	UL16-binding protein 3 (ALCAN-gamma) (NKG2D ligand 3) (N2DL-3) (NKG2DL3) (Retinoic acid early transcript 1N)
--------------------------	--

Note

For Research Use Only , Not for Diagnostic Use

Product Name: Recombinant Human ULBP3 Protein
Catalog #: PHH6508

