



Product Name: Recombinant Human PLGLB1 Protein
Catalog #: PHH6521



Summary

Name	PLGB_HUMAN
Purity	Greater than 85% as determined by reducing SDS-PAGE
Endotoxin level	Please contact with the lab for this information
Tag	N-MBP & C-His
Accession #	Q02325
Host	HEK293
Species	Human
Molecular Mass	51 kDa
Buffer	Phosphate buffered saline
Shipping	The product is shipped on dry ice/polar packs. Upon receipt, store it immediately at the temperature listed below.
Stability&Storage	Store at $\leq -70^{\circ}\text{C}$, stable for 6 months after receipt. Store at $\leq -70^{\circ}\text{C}$, stable for 3 months under sterile conditions after opening. Please minimize freeze-thaw cycles.
Sequence	EPLDDYVNTQ GPSLFSVTKK QLGAGSREEC AAKCEEDKEF TCRAFQYHSK EQQCVIMAEN RKSSIIIRMR DAVLF EK

image



Background

Alternative Names Plasminogen-like protein B (Plasminogen-related protein B)

Note

For Research Use Only , Not for Diagnostic Use