



Product Name: Recombinant Human VAMP5 Protein
Catalog #: PHH6523



Summary

Name	VAMP5_HUMAN
Purity	Greater than 85% as determined by reducing SDS-PAGE
Endotoxin level	Please contact with the lab for this information
Tag	N-MBP & C-His
Accession #	O95183
Host	HEK293
Species	Human
Molecular Mass	51 kDa
Buffer	Phosphate buffered saline
Shipping	The product is shipped on dry ice/polar packs. Upon receipt, store it immediately at the temperature listed below.
Stability&Storage	Store at $\leq -70^{\circ}\text{C}$, stable for 6 months after receipt. Store at $\leq -70^{\circ}\text{C}$, stable for 3 months under sterile conditions after opening. Please minimize freeze-thaw cycles.
Sequence	QQANEVTEIM RNNFGKVLER GVKLAELQQR SDQLDMSST FNKTTQNLAQ KKCWENIRYR

image



Background

Alternative Names Vesicle-associated membrane protein 5 (VAMP-5) (Myobrevin)

Note

For Research Use Only , Not for Diagnostic Use