


Product Name: Recombinant Human UCN2 Protein
Catalog #: PHH6526



Summary

Name	UCN2_HUMAN
Purity	Greater than 85% as determined by reducing SDS-PAGE
Endotoxin level	Please contact with the lab for this information
Tag	C-hFc & C-His
Accession #	Q96RP3
Host	HEK293
Species	Human
Molecular Mass	38 kDa
Buffer	Phosphate buffered saline
Shipping	The product is shipped on dry ice/polar packs. Upon receipt, store it immediately at the temperature listed below.
Stability&Storage	Store at $\leq -70^{\circ}\text{C}$, stable for 6 months after receipt. Store at $\leq -70^{\circ}\text{C}$, stable for 3 months under sterile conditions after opening. Please minimize freeze-thaw cycles.
Sequence	TIPI TFQLRP QNSPQTTPRP AAESPSAAP TWPWAAQSHC SPTRHPGSRI VLSLDVPIGL LQILLEQARA RAAREQATTN ARILARVGHC

image



Background

Alternative Names Urocortin-2 (Stresscopin-related peptide) (Urocortin II) (Ucn II) (Urocortin-related peptide)

Note

For Research Use Only , Not for Diagnostic Use