


Product Name: Recombinant Human ZDHHC9 Protein
Catalog #: PHH6530



Summary

Name	ZDHC9_HUMAN
Purity	Greater than 85% as determined by reducing SDS-PAGE
Endotoxin level	Please contact with the lab for this information
Tag	C-His
Accession #	Q9Y397
Host	HEK293
Species	Human
Molecular Mass	13 kDa
Buffer	Phosphate buffered saline
Shipping	The product is shipped on dry ice/polar packs. Upon receipt, store it immediately at the temperature listed below.
Stability&Storage	Store at $\leq -70^{\circ}\text{C}$, stable for 6 months after receipt. Store at $\leq -70^{\circ}\text{C}$, stable for 3 months under sterile conditions after opening. Please minimize freeze-thaw cycles.
Sequence	ALNQTTNEDI KGSWTGKNRV QNPYSHGNIV KNCCEVLCPG LPPSVLDRRG ILPLEESGSR PPSTQETSSS LLPQSPAPTE HLNSNEMPED SSTPEEMPPP EPPEPPQEAA

image



Background

Alternative Names	Palmitoyltransferase ZDHHC9 (EC 2.3.1.225) (Zinc finger DHHC domain-containing protein 9) (DHHC-9) (DHHC9) (Zinc finger protein 379) (Zinc finger protein 380)
--------------------------	--

Note

For Research Use Only , Not for Diagnostic Use

Product Name: Recombinant Human ZDHHC9 Protein
Catalog #: PHH6530

