


Product Name: Recombinant Human CD226 Protein
Catalog #: PHH6565



Summary

Name	CD226_HUMAN
Purity	Greater than 85% as determined by reducing SDS-PAGE
Endotoxin level	Please contact with the lab for this information
Tag	N-MBP & C-His
Accession #	Q15762
Host	HEK293
Species	Human
Molecular Mass	56 kDa
Buffer	Phosphate buffered saline
Shipping	The product is shipped on dry ice/polar packs. Upon receipt, store it immediately at the temperature listed below.
Stability&Storage	Store at $\leq -70^{\circ}\text{C}$, stable for 6 months after receipt. Store at $\leq -70^{\circ}\text{C}$, stable for 3 months under sterile conditions after opening. Please minimize freeze-thaw cycles.
Sequence	FKIGTQQDSI AIFSPTHGMV IRKPYAERVY FLNSTMASNN MTLFFRNASE DDVGYYSCSL YTPQGTWQK VIQVVQSDSF EAAVPSNSHI VSEPGKNVTL TCQPQMTWPV QAVRWELIQP

image



Background

Alternative Names CD226 antigen (DNAX accessory molecule 1) (DNAM-1) (CD antigen CD226)

Note

For Research Use Only , Not for Diagnostic Use