


Product Name: Recombinant Human ATP12A Protein
Catalog #: PHH6575



Summary

Name	AT12A_HUMAN
Purity	Greater than 85% as determined by reducing SDS-PAGE
Endotoxin level	Please contact with the lab for this information
Tag	C-His
Accession #	P54707
Host	HEK293
Species	Human
Molecular Mass	13 kDa
Buffer	Phosphate buffered saline
Shipping	The product is shipped on dry ice/polar packs. Upon receipt, store it immediately at the temperature listed below.
Stability&Storage	Store at $\leq -70^{\circ}\text{C}$, stable for 6 months after receipt. Store at $\leq -70^{\circ}\text{C}$, stable for 3 months under sterile conditions after opening. Please minimize freeze-thaw cycles.
Sequence	VELSGTKDIV KTDKGDGKEK YRGLKNNCLE LKKKNHKEEF QKELHLDHKL LSNRELEEKY GTDIIMGLSS TRAAELLARD GPNSLTTPPKQ

image



Background

Alternative Names	Potassium-transporting ATPase alpha chain 2 (EC 7.2.2.19) (Non-gastric H(+)/K(+) ATPase subunit alpha) (Proton pump)
--------------------------	--

Note

For Research Use Only , Not for Diagnostic Use