



Product Name: Recombinant Human DEFB119 Protein
Catalog #: PHH6591



Summary

Name	DB119_HUMAN
Purity	Greater than 85% as determined by reducing SDS-PAGE
Endotoxin level	Please contact with the lab for this information
Tag	N-MBP & C-His
Accession #	Q8N690
Host	HEK293
Species	Human
Molecular Mass	50 kDa
Buffer	Phosphate buffered saline
Shipping	The product is shipped on dry ice/polar packs. Upon receipt, store it immediately at the temperature listed below.
Stability&Storage	Store at $\leq -70^{\circ}\text{C}$, stable for 6 months after receipt. Store at $\leq -70^{\circ}\text{C}$, stable for 3 months under sterile conditions after opening. Please minimize freeze-thaw cycles.
Sequence	NSGICRASCK KNEQPYLYCR NCQSCCLQSY MRISISGKEE NTDWSYEKQW

image



Background

Alternative Names	Beta-defensin 119 (Beta-defensin 120) (Beta-defensin 19) (DEFB-19) (Beta-defensin 20) (DEFB-20) (Defensin, beta 119) (Defensin, beta 120) (ESC42-RELA)
--------------------------	--

Note

For Research Use Only , Not for Diagnostic Use