


Product Name: Recombinant Human IL17D Protein
Catalog #: PHH6597



Summary

Name	IL17D_HUMAN
Purity	Greater than 85% as determined by reducing SDS-PAGE
Endotoxin level	Please contact with the lab for this information
Tag	N-MBP & C-His
Accession #	Q8TAD2
Host	HEK293
Species	Human
Molecular Mass	56 kDa
Buffer	Phosphate buffered saline
Shipping	The product is shipped on dry ice/polar packs. Upon receipt, store it immediately at the temperature listed below.
Stability&Storage	Store at $\leq -70^{\circ}\text{C}$, stable for 6 months after receipt. Store at $\leq -70^{\circ}\text{C}$, stable for 3 months under sterile conditions after opening. Please minimize freeze-thaw cycles.
Sequence	AGAPRAGRPP ARPRGCADRP EELLEQLYGR LAAGVLSAFH HTLQLGPREQ ARNASCPAGG RPADRRFRPP TNLRSVSPWA YRISYDPARY PRYLPEAYCL CRGCLTGLFG EEDVRFRSAP VYMPTVVLRR TPACAGGRSV YTEAYVTIPV GCTCVPEPEK DADSINSSID KQGAKLLGP NDAPAGP

image



Background

Alternative Names Interleukin-17D (IL-17D) (Interleukin-27) (IL-27)

Note

For Research Use Only , Not for Diagnostic Use

Product Name: Recombinant Human IL17D Protein
Catalog #: PHH6597

