


Product Name: Recombinant Human IL25 Protein
Catalog #: PHH6609



Summary

Name	IL25_HUMAN
Purity	Greater than 85% as determined by reducing SDS-PAGE
Endotoxin level	Please contact with the lab for this information
Tag	C-His
Accession #	Q9H293
Host	HEK293
Species	Human
Molecular Mass	18 kDa
Buffer	Phosphate buffered saline
Shipping	The product is shipped on dry ice/polar packs. Upon receipt, store it immediately at the temperature listed below.
Stability&Storage	Store at $\leq -70^{\circ}\text{C}$, stable for 12 months after receipt. Store at $\leq -70^{\circ}\text{C}$, stable for 3 months under sterile conditions after opening. Please minimize freeze-thaw cycles.
Sequence	YSHWPSCCPS KGQDTSEELL RWSTVPVPPL EPARPNRHPE SCRASEDGPL NSRAISPWRY ELDRDLNRLP QDLYHARCLC PHCVSLQTGS HMDPRGNSEL LYHNQTVFYR RPCHGEKGTH KGYCLERRLY RVSLACVCVR PRVMG

image



Background

Alternative Names Interleukin-25 (IL-25) (Interleukin-17E) (IL-17E)

Note

For Research Use Only , Not for Diagnostic Use