


Product Name: Recombinant Human IGLON5 Protein
Catalog #: PHH6619



Summary

Name	IGLO5_HUMAN
Purity	Greater than 85% as determined by reducing SDS-PAGE
Endotoxin level	Please contact with the lab for this information
Tag	C-His
Accession #	A6NGN9
Host	HEK293
Species	Human
Molecular Mass	17 kDa
Buffer	Phosphate buffered saline
Shipping	The product is shipped on dry ice/polar packs. Upon receipt, store it immediately at the temperature listed below.
Stability&Storage	Store at $\leq -70^{\circ}\text{C}$, stable for 6 months after receipt. Store at $\leq -70^{\circ}\text{C}$, stable for 3 months under sterile conditions after opening. Please minimize freeze-thaw cycles.
Sequence	EPTVTWRQLR DGFTSEGEIL EISDIQRGQA GEYECVTHNG VNSAPDSRRV LVTVNYPPTI TDVTSARTAL GRAALLRCEA MAVPPADFQW YKDDRLLSSG TAEGLKVQTE RTRSMLLFAN VSARHYGNYT CRAANRLGAS

image



Background

Alternative Names IgLON family member 5

Note

For Research Use Only , Not for Diagnostic Use