



Product Name: Recombinant Human NEGR1 Protein
Catalog #: PHH6625



Summary

Name	NEGR1_HUMAN
Purity	Greater than 85% as determined by reducing SDS-PAGE
Endotoxin level	Please contact with the lab for this information
Tag	C-His
Accession #	Q7Z3B1
Host	HEK293
Species	Human
Molecular Mass	32 kDa
Buffer	Lyophilized from Phosphate buffered saline
Shipping	The product is shipped on dry ice/polar packs. Upon receipt, store it immediately at the temperature listed below.
Stability&Storage	Store at $\leq -70^{\circ}\text{C}$, stable for 12 months after receipt. Store at $\leq -70^{\circ}\text{C}$, stable for 3 months under sterile conditions after opening. Please minimize freeze-thaw cycles.
Sequence	VDFPWAAVDN MMVRKGD TAV LRCYLEDGAS KGAWLNRSSI IFAGGDKWSV DPRVSISTLN KRDYSLQIQN VDVTDDGPYT CSVQTQHTPR TMQVHLTVQV PPKIYDISND MTVNEGTNVT LTCLATGKPE PSISWRHISP SAKPFENGQY LDIYGITRDQ AGEYEC SAEN DVSFPDVRKV KVVVNFAPTI QEIKSGTVTP GRSGLIRCEG AGVPPPAFEW YKGEKKLFNG QQGIIQNFS TRSILVTNV TQEHFGNYTC VAANKLGTTN ASLPLNPPST AQYGITG

image



Background

Alternative Names Neuronal growth regulator 1 (IgLON family member 4)

Product Name: Recombinant Human NEGR1 Protein
Catalog #: PHH6625



Note

For Research Use Only , Not for Diagnostic Use