



Product Name: Recombinant Human DEFB103A Protein
Catalog #: PHH6629



Summary

Name	D103A_HUMAN
Purity	Greater than 85% as determined by reducing SDS-PAGE
Endotoxin level	Please contact with the lab for this information
Tag	N-MBP & C-His
Accession #	P81534
Host	HEK293
Species	Human
Molecular Mass	48 kDa
Buffer	Phosphate buffered saline
Shipping	The product is shipped on dry ice/polar packs. Upon receipt, store it immediately at the temperature listed below.
Stability&Storage	Store at $\leq -70^{\circ}\text{C}$, stable for 6 months after receipt. Store at $\leq -70^{\circ}\text{C}$, stable for 3 months under sterile conditions after opening. Please minimize freeze-thaw cycles.
Sequence	YYCRVRGGRC AVLSCLPKEE QIGKCSTRGR

image



Background

Alternative Names Beta-defensin 103 (Beta-defensin 3) (BD-3) (DEFB-3) (HBD3) (hBD-3) (Defensin, beta 103) (Defensin-like protein)

Note

For Research Use Only , Not for Diagnostic Use