


Product Name: Recombinant Human NAALADL1 Protein
Catalog #: PHH6642



Summary

| | |
|------------------------------|--|
| Name | NALDL_HUMAN |
| Purity | Greater than 85% as determined by reducing SDS-PAGE |
| Endotoxin level | Please contact with the lab for this information |
| Tag | C-His |
| Accession # | Q9UQQ1 |
| Host | HEK293 |
| Species | Human |
| Molecular Mass | 12 kDa |
| Buffer | Phosphate buffered saline |
| Shipping | The product is shipped on dry ice/polar packs. Upon receipt, store it immediately at the temperature listed below. |
| Stability&Storage | Store at $\leq -70^{\circ}\text{C}$, stable for 6 months after receipt. Store at $\leq -70^{\circ}\text{C}$, stable for 3 months under sterile conditions after opening. Please minimize freeze-thaw cycles. |
| Sequence | DLEILETVMG QLDAHRIREN LRELSREPHL ASSPRDEDLV QLLLQRWKDP ESGLDSAEAS TYEVLLSFPS QEQPNVVDIV GPTGGIIHSC HRTEENVGTGE |

image



Background

| | |
|--------------------------|--|
| Alternative Names | Aminopeptidase NAALADL1 (EC 3.4.11.-) (100 kDa ileum brush border membrane protein) (I100) (Ileal dipeptidylpeptidase) (N-acetylated-alpha-linked acidic dipeptidase-like protein) (NAALADase L) |
|--------------------------|--|

Note

For Research Use Only , Not for Diagnostic Use