


Product Name: Recombinant Human CD160 Protein
Catalog #: PHH6655



Summary

Name	BY55_HUMAN
Purity	Greater than 85% as determined by reducing SDS-PAGE
Endotoxin level	Please contact with the lab for this information
Tag	N-His & N-MBP
Accession #	O95971
Host	HEK293
Species	Human
Molecular Mass	57 kDa
Buffer	Phosphate buffered saline
Shipping	The product is shipped on dry ice/polar packs. Upon receipt, store it immediately at the temperature listed below.
Stability&Storage	Store at $\leq -70^{\circ}\text{C}$, stable for 6 months after receipt. Store at $\leq -70^{\circ}\text{C}$, stable for 3 months under sterile conditions after opening. Please minimize freeze-thaw cycles.
Sequence	SSASQEGTRL NLICTVWHKK EEAEGFVVFL CKDRSGDCSP ETSLKQLRLK RDPGIDGVGE ISSQLMFTIS QVTPLHSGTY QCCARSQKSG IRLQGHFFSI LFTETGNYTV TGLKQRQHLE

image



Background

Alternative Names CD160 antigen (Natural killer cell receptor BY55) (CD antigen CD160)
[Cleaved into: CD160 antigen, soluble form]

Note

For Research Use Only , Not for Diagnostic Use