


Product Name: Recombinant Human EFNA4 Protein
Catalog #: PHH6667



Summary

Name	EFNA4_HUMAN
Purity	Greater than 85% as determined by reducing SDS-PAGE
Endotoxin level	Please contact with the lab for this information
Tag	C-hFc & C-His
Accession #	P52798
Host	HEK293
Species	Human
Molecular Mass	44 kDa
Buffer	Phosphate buffered saline
Shipping	The product is shipped on dry ice/polar packs. Upon receipt, store it immediately at the temperature listed below.
Stability&Storage	Store at $\leq -70^{\circ}\text{C}$, stable for 6 months after receipt. Store at $\leq -70^{\circ}\text{C}$, stable for 3 months under sterile conditions after opening. Please minimize freeze-thaw cycles.
Sequence	LRHVYWNSS NPRLRGDAV VELGLNDYLD IVCPHYEGPG PPEGPETFAL YMVDWPGYES CQAEGRPRAYK RWVCSLPFGH VQFSEKIQRFP TFSLGFEFL PGETYYYISV PTPESGQCL RLQVSVCCKE RKSESAHPVG SPGESG

image



Background

Alternative Names Ephrin-A4 (EPH-related receptor tyrosine kinase ligand 4) (LERK-4)

Note

For Research Use Only , Not for Diagnostic Use