

**Product Name: Recombinant Human ISG15 Protein**  
**Catalog #: PHH6669**



## Summary

<b>Name</b>	ISG15_HUMAN
<b>Purity</b>	Greater than 85% as determined by reducing SDS-PAGE
<b>Endotoxin level</b>	Please contact with the lab for this information
<b>Tag</b>	C-His
<b>Accession #</b>	P05161
<b>Host</b>	HEK293
<b>Species</b>	Human
<b>Molecular Mass</b>	18 kDa
<b>Buffer</b>	Phosphate buffered saline
<b>Shipping</b>	The product is shipped on dry ice/polar packs. Upon receipt, store it immediately at the temperature listed below.
<b>Stability&amp;Storage</b>	Store at $\leq -70^{\circ}\text{C}$ , stable for 6 months after receipt. Store at $\leq -70^{\circ}\text{C}$ , stable for 3 months under sterile conditions after opening. Please minimize freeze-thaw cycles.
<b>Sequence</b>	MGWDLTVKML AGNEFQVSLSSSMSVSELKA QITQKIGVHA FQQLAVHPS GVALQDRVPL ASQGLGPGST VLLVVDKCDE PLSILVRNNGRSTYEVRL TQTV AHLKQQ VSGLEGVQDD LFWLTFEGKP LEDQLPLGEY GLKPLSTVFM NLRLRGG

## image



## Background

**Alternative Names** Ubiquitin-like protein ISG15 (Interferon-induced 15 kDa protein) (Interferon-induced 17 kDa protein) (IP17) (Ubiquitin cross-reactive protein) (hUCRP)

## Note

For Research Use Only , Not for Diagnostic Use

**Product Name: Recombinant Human ISG15 Protein**  
**Catalog #: PHH6669**

---

