



Product Name: Recombinant Human ADM Protein
Catalog #: PHH6680



Summary

Name	ADML_HUMAN
Purity	Greater than 85% as determined by reducing SDS-PAGE
Endotoxin level	Please contact with the lab for this information
Tag	N-His & N-hFc
Accession #	P35318
Host	HEK293
Species	Human
Molecular Mass	34 kDa
Buffer	Phosphate buffered saline
Shipping	The product is shipped on dry ice/polar packs. Upon receipt, store it immediately at the temperature listed below.
Stability&Storage	Store at $\leq -70^{\circ}\text{C}$, stable for 6 months after receipt. Store at $\leq -70^{\circ}\text{C}$, stable for 3 months under sterile conditions after opening. Please minimize freeze-thaw cycles.
Sequence	NFQGLRSFGC RFGTCTVQKL AHQIQFTDK DKDNVAPRSK

image



Background

Alternative Names Pro-adrenomedullin [Cleaved into: Adrenomedullin (AM); Proadrenomedullin N-20 terminal peptide (ProAM N-terminal 20 peptide) (PAMP) (ProAM-N20)]

Note

For Research Use Only , Not for Diagnostic Use