


Product Name: Recombinant Human CALCR Protein
Catalog #: PHH6683



Summary

Name	CALCR_HUMAN
Purity	Greater than 85% as determined by reducing SDS-PAGE
Endotoxin level	Please contact with the lab for this information
Tag	C-hFc & C-His
Accession #	P30988
Host	HEK293
Species	Human
Molecular Mass	43 kDa
Buffer	Phosphate buffered saline
Shipping	The product is shipped on dry ice/polar packs. Upon receipt, store it immediately at the temperature listed below.
Stability&Storage	Store at $\leq -70^{\circ}\text{C}$, stable for 6 months after receipt. Store at $\leq -70^{\circ}\text{C}$, stable for 3 months under sterile conditions after opening. Please minimize freeze-thaw cycles.
Sequence	YPTIEPKPFL YVVGRRKKMMD AQYKCYDRMQ QLPAYQGEGP YCNRTWDGWL CWDDTPAGVL SYQFCPDYFP DFDPSEKVTK YCDEKGVWFK HPENNRTWSN YTMCNFTPE KLKNAYVLYY

image



SDS



WB

Background

Alternative Names Calcitonin receptor (CT-R)

Note

For Research Use Only , Not for Diagnostic Use