



Product Name: Recombinant Human SLC1A7 Protein
Catalog #: PHH6688



Summary

Name	EAA5_HUMAN
Purity	Greater than 85% as determined by reducing SDS-PAGE
Endotoxin level	Please contact with the lab for this information
Tag	C-hFc & C-His
Accession #	O00341
Host	HEK293
Species	Human
Molecular Mass	39 kDa
Buffer	Phosphate buffered saline
Shipping	The product is shipped on dry ice/polar packs. Upon receipt, store it immediately at the temperature listed below.
Stability&Storage	Store at $\leq -70^{\circ}\text{C}$, stable for 12 months after receipt. Store at $\leq -70^{\circ}\text{C}$, stable for 3 months under sterile conditions after opening. Please minimize freeze-thaw cycles.
Sequence	QKETTEQSGK PIMSSADALL DLIRNMFAN LVEATFKQYR TKTTPVVKSP KVAPEEAPPR RILYGVQEE NGSHVQNFAL DLTPPPEVVY

image



Background

Alternative Names Excitatory amino acid transporter 5 (Retinal glutamate transporter) (Solute carrier family 1 member 7)

Note

For Research Use Only , Not for Diagnostic Use