


Summary

Name	ECE1_HUMAN
Purity	Greater than 85% as determined by reducing SDS-PAGE
Endotoxin level	Please contact with the lab for this information
Tag	N-9×His
Accession #	P42892
Host	HEK293
Species	Human
Molecular Mass	79 kDa
Buffer	Phosphate buffered saline
Shipping	The product is shipped on dry ice/polar packs. Upon receipt, store it immediately at the temperature listed below.
Stability&Storage	Store at $\leq -70^{\circ}\text{C}$, stable for 6 months after receipt. Store at $\leq -70^{\circ}\text{C}$, stable for 3 months under sterile conditions after opening. Please minimize freeze-thaw cycles.
Sequence	YQTRSPSVCL SEACVSVTSS ILSSMDPTVD PCHDFFSYAC GGWIKANPVP DGHSRWGTFS NLWEHNQAI KHLLNSTAS VSEAERKAQV YRACMNETR IEELRAKPLM ELIERLGGWN ITGPWAKDNF QDTLQVVTAH YRTSPFFSVY VSADSKNSNS NVIQVDQSGL GLPSRDYYLN KTENEKVLTG YLNYMVQLGK LLGGGDEEAI RPQMQQILDF ETALANITIP QEKRREDELI YHKVTAELQ TLAPAINWLP FLNTIFYPVE INESEPIVVY DKEYLEQIST LINTTDRCLL NNYMIWNLVR KTSSFLDQRF QDADEKFMEV MYGTTKTCPL RWKFCVSDTE NNLGFALGPM FVKATFAEDS KSIATEIILE IKKAFFESLS TLKWMDEETR KSAKEKADAI YNMIGYPNFI MDPKELDKVF NDYTAVPDLY FENAMRFFNF SWRVTADQLR KAPNRDQWSM TPPMVNAYYS PTKNEIVFPA GILQAPFYTR SSPKALNFGG IGVVVGHELT HAFDDQGREY DKDGNLRPWW KNSSVEAFKR QTECMVEQYS NYSVNGEPVN GRHTLGENIA DNGGLKAAYR AYQNWVKKNG AEHSLPTLGL TNNQLFFLGF AQVWCSVRTP ESSHEGLITD PHSPSRFRVI GSLSNSKEFS EHFRCPPGSP MNPPHKCEVW

image



Product Name: Recombinant Human ECE1 Protein
Catalog #: PHH6689



SDS

WB

Background

Alternative Names

Endothelin-converting enzyme 1 (ECE-1) (EC 3.4.24.71)

Note

For Research Use Only , Not for Diagnostic Use