



Product Name: Recombinant Human BTNL9 Protein
Catalog #: PHH6702



Summary

Name	BTNL9_HUMAN
Purity	Greater than 85% as determined by reducing SDS-PAGE
Endotoxin level	Please contact with the lab for this information
Tag	C-hFc & C-His
Accession #	Q6UXG8
Host	HEK293
Species	Human
Molecular Mass	53 kDa
Buffer	Phosphate buffered saline
Shipping	The product is shipped on dry ice/polar packs. Upon receipt, store it immediately at the temperature listed below.
Stability&Storage	Store at $\leq -70^{\circ}\text{C}$, stable for 6 months after receipt. Store at $\leq -70^{\circ}\text{C}$, stable for 3 months under sterile conditions after opening. Please minimize freeze-thaw cycles.
Sequence	PGEPSSEVKV LGPEYPILAL VGEEVEFPCH LWPQLDAQQM EIRWFRSQTF NVVHLYQEQQ ELPGRQMPAF RNRTKLVKDD IAYGSVVLQL HSIIPSDKGT YGRFHSDFN SGEALWELEV AGLGSDPHLS LEGFKEGGIQ LRLRSSGWYP KPKVQWRDHQ GQCLPPEFEA IVWDAQDLFS LETSVVVRAG ALSNVSVSIQ NLLSQKKEK VVQIADVFP

image



Background

Alternative Names Butyrophilin-like protein 9

Note

Product Name: Recombinant Human BTNL9 Protein
Catalog #: PHH6702



For Research Use Only , Not for Diagnostic Use