



Product Name: Recombinant Human CLEC12A Protein
Catalog #: PHH6716



Summary

Name	CL12A_HUMAN
Purity	Greater than 85% as determined by reducing SDS-PAGE
Endotoxin level	Please contact with the lab for this information
Tag	N-MBP & C-His
Accession #	Q5QGZ9
Host	HEK293
Species	Human
Molecular Mass	65 kDa
Buffer	Phosphate buffered saline
Shipping	The product is shipped on dry ice/polar packs. Upon receipt, store it immediately at the temperature listed below.
Stability&Storage	Store at $\leq -70^{\circ}\text{C}$, stable for 6 months after receipt. Store at $\leq -70^{\circ}\text{C}$, stable for 3 months under sterile conditions after opening. Please minimize freeze-thaw cycles.
Sequence	ISEELQRNIS LQLMSNMNIS NKIRNLSTTL QTIATKLCRE LYSKEQEHC KPCPRRWIWH KDSCYFLSDD VQTWQESKMA CAAQNASLLK INNKNALEFI KSQSRSDYDW LGLSPEEDST RGMVRVDNIIN SSAWVIRNAP DLNNMYCGYI NRLYVQYYHC TYKKRMICEK MANPVQLGST

image



Background

Alternative Names	C-type lectin domain family 12 member A (C-type lectin-like molecule 1) (CLL-1) (Dendritic cell-associated lectin 2) (DCAL-2) (Myeloid inhibitory C-type lectin-like receptor) (MICL) (CD antigen CD371)
--------------------------	--

Note

Product Name: Recombinant Human CLEC12A Protein
Catalog #: PHH6716



For Research Use Only , Not for Diagnostic Use