



Product Name: Recombinant Human GHRL Protein
Catalog #: PHH6740



Summary

Name	GHRL_HUMAN
Purity	Greater than 85% as determined by reducing SDS-PAGE
Endotoxin level	Please contact with the lab for this information
Tag	N-9×His
Accession #	Q9UBU3
Host	HEK293
Species	Human
Molecular Mass	12 kDa
Buffer	Phosphate buffered saline
Shipping	The product is shipped on dry ice/polar packs. Upon receipt, store it immediately at the temperature listed below.
Stability&Storage	Store at $\leq -70^{\circ}\text{C}$, stable for 6 months after receipt. Store at $\leq -70^{\circ}\text{C}$, stable for 3 months under sterile conditions after opening. Please minimize freeze-thaw cycles.
Sequence	EHQRVQQRKE SKKPPAKLQP RALAGWLRPE DGGQAE GAED ELEVRFNAPF DVGIKLSGVQ YQQHSQALGK FLQDILWEEA

image



Background

Alternative Names	Appetite-regulating hormone (Growth hormone secretagogue) (Growth hormone-releasing peptide) (Motilin-related peptide) (Protein M46) [Cleaved into: Ghrelin-27; Ghrelin-28 (Ghrelin); Obestatin]
--------------------------	--

Note

For Research Use Only , Not for Diagnostic Use