



Product Name: Recombinant Human IFI30 Protein
Catalog #: PHH6743



Summary

Name	GILT_HUMAN
Purity	Greater than 85% as determined by reducing SDS-PAGE
Endotoxin level	Please contact with the lab for this information
Tag	N-9×His
Accession #	P13284
Host	HEK293
Species	Human
Molecular Mass	27 kDa
Buffer	Phosphate buffered saline
Shipping	The product is shipped on dry ice/polar packs. Upon receipt, store it immediately at the temperature listed below.
Stability&Storage	Store at $\leq -70^{\circ}\text{C}$, stable for 6 months after receipt. Store at $\leq -70^{\circ}\text{C}$, stable for 3 months under sterile conditions after opening. Please minimize freeze-thaw cycles.
Sequence	SPLQALDFFG NGPPVNYKTG NLYLRGPLKK SNAPLVNVTL YYEALCGGCR AFLIRELFPT WLLVMEILNV TLVPYGNAQE QNVSGRWEFK CQHGEEECKF NKVEACVLDE LDMELAFIT VCMEEFEDME RSLPLCLQLY APGLSPDTIM ECAMGDRGMQ LMHANAQRTD ALQPPHEYVP WVTVNGKPLE DQTQLLTLVC QLYQGKKPDV CPSSTSSLRS VCFK

image



Background

Alternative Names	Gamma-interferon-inducible lysosomal thiol reductase (EC 1.8.-.-) (Gamma-interferon-inducible protein IP-30) (Legumaturain)
--------------------------	---

Note

Product Name: Recombinant Human IFI30 Protein
Catalog #: PHH6743



For Research Use Only , Not for Diagnostic Use