


Product Name: Recombinant Human ADAM19 Protein
Catalog #: PHH6750



Summary

Name	ADA19_HUMAN
Purity	Greater than 85% as determined by reducing SDS-PAGE
Endotoxin level	Please contact with the lab for this information
Tag	C-His
Accession #	Q9H013
Host	HEK293
Species	Human
Molecular Mass	10 kDa
Buffer	Phosphate buffered saline
Shipping	The product is shipped on dry ice/polar packs. Upon receipt, store it immediately at the temperature listed below.
Stability&Storage	Store at $\leq -70^{\circ}\text{C}$, stable for 6 months after receipt. Store at $\leq -70^{\circ}\text{C}$, stable for 3 months under sterile conditions after opening. Please minimize freeze-thaw cycles.
Sequence	FEKVNvagdt FGNCgKDMNG EHRKCNMRDA KCGKIQCQSS EARPLESNAV PIDTTIIMNG RQIQCRGTHV

image



Background

Alternative Names	Disintegrin and metalloproteinase domain-containing protein 19 (ADAM 19) (EC 3.4.24.-) (Meltrin-beta) (Metalloprotease and disintegrin dendritic antigen marker) (MADDAM)
--------------------------	---

Note

For Research Use Only , Not for Diagnostic Use